


## PPAR $\gamma$ Ligand Screening/Characterization Assay Kit (Fluorometric) (#BN00680)

(Catalog # BN00680, 100 assays; Store kit at -20°C)

### I. Introduction:

The Peroxisome Proliferator Activated Receptor (PPAR) family of ligand-activated transcription factors consists of three subtypes encoded by separate genes: PPAR $\alpha$ , PPAR $\delta$  and PPAR $\gamma$ . Of these, PPAR $\gamma$  plays an important role in the regulation of fatty acid storage and glucose metabolism. The genes activated by PPAR $\gamma$  stimulate lipid uptake and adipogenesis by fat cells. Many endogenous molecules such as, polyunsaturated fatty acids like arachidonic acid and its metabolites, are known to bind and activate PPAR $\gamma$ . The binding of activating ligands to the ligand binding domain (LBD) of PPAR $\gamma$  promotes its heterodimerization with retinoic acid-like receptor (RXR), which results in the regulated expression of target genes involved in lipid metabolism. Such ligand-based activation of PPAR $\gamma$  may be responsible for inhibiting the growth of cultured human breast, gastric, lung, prostate and other cancer cell lines. In addition, the thiazolidinedione-based anti-diabetic drugs activate PPAR $\gamma$  with greater specificity than PPAR $\alpha$ . Assay Genie's PPAR $\gamma$  Ligand Screening Assay Kit provides a single step fluorescence-based assay for screening potential PPAR $\gamma$ -specific ligands. The assay utilizes the ability of potential PPAR $\gamma$ -binding ligands to displace a fluorescent probe, which has a strong affinity for PPAR $\gamma$  Ligand Binding Domain, resulting in loss of fluorescence of the probe. The relative drop in the fluorescence, as a result of competitive binding of PPAR $\gamma$  ligand, can be correlated to the affinity (and hence IC<sub>50</sub>) of the PPAR $\gamma$  candidate ligand. Assay Genie's PPAR $\gamma$  Ligand Screening Assay Kit is easy to use, faster and more convenient as compared to Fluorescence Polarization and TR-FRET-based screening methods. The assay kit can be used to identify and characterize PPAR $\gamma$ -specific ligands for therapeutic applications.



### II. Applications:

- Screening of potential PPAR $\gamma$  binding ligands.

### III. Kit Contents:

Components	BN00680	Cap Code	Part Number
PPAR $\gamma$ Assay Buffer	25 ml	WM	BN00680-1
PPAR $\gamma$ Assay Probe	10 $\mu$ l	Red	BN00680-2
PPAR $\gamma$ (Human Recombinant)	2 x 250 $\mu$ l	Brown	BN00680-3
PPAR $\gamma$ Ligand Control (100 mM in DMSO)	10 $\mu$ l	Blue	BN00680-4
384-well Low Volume Black Plate	1 Plate	-	BN00680-5

### IV. User Supplied Reagents and Equipment:

- DMSO, 384-well black plate.
- Multi-well spectrofluorometer.

### V. Storage Conditions and Reagent Preparation:

Store kit at -20°C, protected from light. Briefly centrifuge small vials at low speed prior to opening. Read the entire protocol before performing the assay.

- PPAR $\gamma$  Assay Buffer:** Bring to room temperature before use. Store at -20°C. Avoid prolonged storage of the PPAR $\gamma$  Assay Buffer at room temperature or 4°C.
- Human PPAR $\gamma$ :** Store at -80°C. Avoid repeated freeze/thaw cycles. Each vial contains enough protein for 50 assays.
- PPAR $\gamma$  Assay Probe and Ligand Control:** Store at -20°C. Bring to room temperature before use.

### VI. PPAR $\gamma$ Ligand Screening Assay Protocol:

- PPAR $\gamma$  Assay probe preparation:** Dilute 5  $\mu$ l of the PPAR $\gamma$  Assay Probe with 495  $\mu$ l of DMSO. Mix well by light Vortexing. Use the probe immediately.
- Screening Compounds, Inhibitor Control & Blank Control Preparations:** Dissolve the test ligands in DMSO or other appropriate solvent. Use 1  $\mu$ l of test ligand (Sample, **S**) or 1  $\mu$ l DMSO (Solvent Control, **SC**) into empty well(s). For Ligand Control (**LC**), dilute 10X by adding 1  $\mu$ l of PPAR $\gamma$  Ligand Control to 9  $\mu$ l DMSO. Use 1  $\mu$ l of 10x diluted PPAR $\gamma$  Ligand Control (in DMSO) into each well(s). In order to obtain IC<sub>50</sub> values, different concentrations of test ligand and/or PPAR $\gamma$  Ligand Control should be tested.
- PPAR $\gamma$  Assay Mix:** Based on number of samples to be tested, prepare appropriate amount of PPAR $\gamma$  Assay Mix per well as below:

PPAR $\gamma$ Protein	5 $\mu$ l
PPAR $\gamma$ Assay Probe (diluted)	1 $\mu$ l
PPAR $\gamma$ Assay Buffer	18 $\mu$ l
Total Volume	24 $\mu$ l

Mix well by pipetting up and down. Incubate at RT for 5-10 min. Add 24  $\mu$ l of PPAR $\gamma$  Assay Mix to each well containing test, solvent control and ligand control. Incubate at RT for 5 min before reading. Final reaction volume in each well shouldn't exceed 25  $\mu$ l. Store unused PPAR $\gamma$  protein immediately at -80°C.


#### Notes:

- If the test ligand is insoluble at high concentrations, precipitation might be observed during the assay. In that case, DMSO can be used up to 10% of final assay volume to increase the solubility of the test ligand in final assay solution.

!

4. **Measurement:** Measure the fluorescence intensity (Ex/Em = 375/460-470 nm) of the samples and the controls in an endpoint mode. The fluorescence signal is stable up to 1 h with minimum loss.
5. **Calculations:** Plot the % Relative Fluorescence (RFU, drop in the fluorescence intensity) and plot it against increasing concentration of the test ligand in the assay as given below; obtain IC<sub>50</sub>.

$$\% \text{ Relative Fluorescence} = \frac{RFU (S)}{RFU (SC)} \times 100$$



**Figure:** A variety of PPAR $\gamma$ -specific ligands (GW 1929, Rosiglitazone, Ciglitazone and Pioglitazone) and a PPAR $\alpha$ -specific ligand (WY 14643) were tested using PPAR $\gamma$  Ligand Screening Assay Kit. Assays were performed following the kit protocol.

**FOR RESEARCH USE ONLY! Not to be used on humans.**