

## Xanthine Oxidase Activity Colorimetric/Fluorometric Assay Kit (#BN00930)

(Catalog #BN00930; 100 assays; Store kit at -20° C)

### I. Introduction:

Xanthine oxidase (XO, EC 1.17.3.2) is present in appreciable amounts in the liver and jejunum in healthy individuals. However, in various liver disorders, XO is released into circulation. Therefore, determination of serum XO level serves as a sensitive indicator of acute liver damage such as jaundice. Assay Genie has developed an easy and sensitive assay to determine XO in variety of samples. In the assay, XO oxidizes xanthine to hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) which reacts stoichiometrically with GenieRed Probe to generate color (at λ = 570 nm) and fluorescence (at Ex/Em = 535/587 nm). Since the color or fluorescence intensity is proportional to XO content, the XO activity can be accurately measured. The kit detects 1-100 mU xanthine oxidase in 100 μl reaction volume.

### II. Kit Contents:

Components	BN00930	Cap Code	Part Number
XO Assay Buffer	25 ml	WM	BN00930-1
GenieRed Probe (in DMSO)	200 μl	Red	BN00930-2A
XO Enzyme Mix	lyophilized	Green	BN00930-4
XO Substrate Mix	lyophilized	Purple	BN00930-5
XO Positive Control	8 μl	Blue	BN00930-6
H <sub>2</sub> O <sub>2</sub> Standard (0.88 M)	0.1 ml	Yellow	BN00930-7

### III. Reagent Preparation and Storage Conditions:

**GenieRed Probe:** Ready to use as supplied. (Need to warm > 20°C to melt frozen DMSO). Store at -20°C, use within two months.

**XO Enzyme Mix:** Dissolve with 220 μl dH<sub>2</sub>O. Pipette up and down to dissolve completely.

**XO Substrate Mix:** Dissolve with 220 μl dH<sub>2</sub>O. Pipette up and down to dissolve completely.

**XO Positive Control:** Dilute with 92 μl dH<sub>2</sub>O. Pipette up and down to dissolve completely. All components in kit should store at -20° C and use within two months.

### IV. Xanthine Oxidase Assay Protocol:

#### 1. Standard Curve Preparations:

Dilute 4 μl of 0.88 M H<sub>2</sub>O<sub>2</sub> Standard into 348 μl dH<sub>2</sub>O to generate 10 mM H<sub>2</sub>O<sub>2</sub> Standard, then dilute 20 μl of 10 mM H<sub>2</sub>O<sub>2</sub> Standard into 980 μl dH<sub>2</sub>O to generate 0.2 mM H<sub>2</sub>O<sub>2</sub> Standard.

**Colorimetric assay:** Add 0, 10, 20, 30, 40, 50 μl of the 0.2 mM H<sub>2</sub>O<sub>2</sub> Standard into 96-well plate in duplicates, bring the total volume to 50 μl each well with dH<sub>2</sub>O to generate 0, 2, 4, 6, 8, 10 nmol/well H<sub>2</sub>O<sub>2</sub> Standard.

**Fluorometric assay:** Dilute 50 μl fresh 0.2 mM H<sub>2</sub>O<sub>2</sub> into 950 μl dH<sub>2</sub>O to generate 10 μM H<sub>2</sub>O<sub>2</sub> Standard. Add 0, 10, 20, 30, 40, 50 μl of the 10 μM H<sub>2</sub>O<sub>2</sub> into 96-well plate in duplicates, bring volume to 50 μl with dH<sub>2</sub>O to generate 0, 0.1, 0.2, 0.3, 0.4, 0.5 nmol/well H<sub>2</sub>O<sub>2</sub> Standard.

**2. Sample and Positive Control Preparations:** Prepare test samples in 50 μl/well with assay buffer in a 96-well plate. Serum can be directly added into sample wells, and adjust volume to 50 μl/well with dH<sub>2</sub>O. Tissues or cells can be extracted with 4 volumes of the Assay Buffer, centrifuge (16,000 x g, 10 min) to get clear XO extract. For the positive control, add 5 μl positive control solution to wells, adjust volume to 50 μl/well with dH<sub>2</sub>O. H<sub>2</sub>O<sub>2</sub> in the sample will generate background. It is important to set up a background control. We suggest using several doses of your sample to ensure the readings are within the linear range.

**3. Reaction Mix Preparation:** Mix enough reagents for the number of assays and standard to be performed. For each well, prepare a total 50 μl Reaction Mix containing:

Xanthine Oxidase Measurement	Background Control
44 μl Assay Buffer	46 μl Assay Buffer
2 μl Substrate Mix	-----
2 μl Enzyme Mix	2 μl Enzyme Mix
2 μl GenieRed Probe**	2 μl GenieRed Probe

\*\* For the fluorescent assay, dilute GenieRed probe 10X to reduce background readings.

**4.** Add 50 μl of the reaction mix to each well containing the H<sub>2</sub>O<sub>2</sub> Standard, Positive Control, and test samples, mix well.

**5.** Measure the plate immediately (OD = 570 nm for colorimetric assay or at Ex/Em = 535/587 nm for fluorometric assay) at T<sub>1</sub> to read A<sub>1</sub>, measure again at T<sub>2</sub> after incubating the reaction at 25° C for 10 - 20 min (or incubate longer time if the sample XO activity is low) to read A<sub>2</sub>, protect from light. The signal generated by XO is ΔA = A<sub>2</sub> - A<sub>1</sub>

#### Notes:

**1)** It is essential to read A<sub>1</sub> and A<sub>2</sub> in the reaction linear range. It will be more accurate if you read the reaction kinetics. Then choose A<sub>1</sub> and A<sub>2</sub> in the reaction linear range.

**2)** Read H<sub>2</sub>O<sub>2</sub> standard after 20 min incubation without subtract A<sub>1</sub>. The standard is stable for a few hours.

**6. Calculation:** Subtract background from all readings. Plot the H<sub>2</sub>O<sub>2</sub> standard Curve. Apply sample ΔA to the H<sub>2</sub>O<sub>2</sub> standard curve to get B nmol of H<sub>2</sub>O<sub>2</sub> (H<sub>2</sub>O<sub>2</sub> generated between T<sub>1</sub> and T<sub>2</sub> in the reaction by XO).

$$\text{XO Activity} = \frac{B}{(T_2 - T_1) \times V} \times \text{Sample Dilution Factor} = \text{nmol/min/ml} = \text{mU/ml}$$


**Where:** B is the amount of H<sub>2</sub>O<sub>2</sub> generated by XO from standard curve (in nmol).

T<sub>1</sub> is the time of the first reading (A<sub>1</sub>) (in min).

T<sub>2</sub> is the time of the second reading (A<sub>2</sub>) (in min).

V is the pretreated sample volume added into the reaction well (in ml).

**Unit Definition:** One unit xanthine oxidase is defined as the amount of enzyme catalyzes the oxidation of xanthine, yielding 1.0 μmol of uric acid and H<sub>2</sub>O<sub>2</sub> per minute at 25° C.



**FOR RESEARCH USE ONLY! Not to be used on humans.**

## GENERAL TROUBLESHOOTING GUIDE:

Problems	Cause	Solution
Assay not working	<ul style="list-style-type: none"> <li>• Use of ice-cold assay buffer</li> <li>• Omission of a step in the protocol</li> <li>• Plate read at incorrect wavelength</li> <li>• Use of a different 96-well plate</li> </ul>	<ul style="list-style-type: none"> <li>• Assay buffer must be at room temperature</li> <li>• Refer and follow the data sheet precisely</li> <li>• Check the wavelength in the data sheet and the filter settings of the instrument</li> <li>• Fluorescence: Black plates (clear bottoms) ; Luminescence: White plates ; Colorimeters: Clear plates</li> </ul>
Samples with erratic readings	<ul style="list-style-type: none"> <li>• Use of an incompatible sample type</li> <li>• Samples prepared in a different buffer</li> <li>• Cell/ tissue samples were not completely homogenized</li> <li>• Samples used after multiple free-thaw cycles</li> <li>• Presence of interfering substance in the sample</li> <li>• Use of old or inappropriately stored samples</li> </ul>	<ul style="list-style-type: none"> <li>• Refer data sheet for details about incompatible samples</li> <li>• Use the assay buffer provided in the kit or refer data sheet for instructions</li> <li>• Use Dounce homogenizer (increase the number of strokes); observe for lysis under microscope</li> <li>• Aliquot and freeze samples if needed to use multiple times</li> <li>• Troubleshoot if needed</li> <li>• Use fresh samples or store at correct temperatures until use</li> </ul>
Lower/ Higher readings in Samples and Standards	<ul style="list-style-type: none"> <li>• Improperly thawed components</li> <li>• Use of expired kit or improperly stored reagents</li> <li>• Allowing the reagents to sit for extended times on ice</li> <li>• Incorrect incubation times or temperatures</li> <li>• Incorrect volumes used</li> </ul>	<ul style="list-style-type: none"> <li>• Thaw all components completely and mix gently before use</li> <li>• Always check the expiry date and store the components appropriately</li> <li>• Always thaw and prepare fresh reaction mix before use</li> <li>• Refer datasheet &amp; verify correct incubation times and temperatures</li> <li>• Use calibrated pipettes and aliquot correctly</li> </ul>
Readings do not follow a linear pattern for Standard curve	<ul style="list-style-type: none"> <li>• Use of partially thawed components</li> <li>• Pipetting errors in the standard</li> <li>• Pipetting errors in the reaction mix</li> <li>• Air bubbles formed in well</li> <li>• Standard stock is at an incorrect concentration</li> <li>• Calculation errors</li> <li>• Substituting reagents from older kits/ lots</li> </ul>	<ul style="list-style-type: none"> <li>• Thaw and resuspend all components before preparing the reaction mix</li> <li>• Avoid pipetting small volumes</li> <li>• Prepare a master reaction mix whenever possible</li> <li>• Pipette gently against the wall of the tubes</li> <li>• Always refer the dilutions in the data sheet</li> <li>• Recheck calculations after referring the data sheet</li> <li>• Use fresh components from the same kit</li> </ul>
Unanticipated results	<ul style="list-style-type: none"> <li>• Measured at incorrect wavelength</li> <li>• Samples contain interfering substances</li> <li>• Use of incompatible sample type</li> <li>• Sample readings above/below the linear range</li> </ul>	<ul style="list-style-type: none"> <li>• Check the equipment and the filter setting</li> <li>• Troubleshoot if it interferes with the kit</li> <li>• Refer data sheet to check if sample is compatible with the kit or optimization is needed</li> <li>• Concentrate/ Dilute sample so as to be in the linear range</li> </ul>
<p><b>Note:</b> The most probable list of causes is under each problem section. Causes/ Solutions may overlap with other problems.</p>		