

O-GlcNAc Modified Glycoprotein Assay Kit (FACS/Microscopy, Green Fluorescence) (#BN00934)
(Catalog # BN00934; 100 assays; Store at -20°C)

I. Introduction:

Glycans are vital components of glycoproteins, glycolipids, and proteoglycans in all domains of life. Glycoproteins are grouped by the type of carbohydrate and amino acid linkage site. *N*-linked glycosylation is a modification of asparagine, whereas *O*-linked glycosylation occurs through the hydroxyl group of serine and threonine residues. Glycosylation occurs co- or post-translationally on >50% of eukaryotic proteins resulting in membrane-associated, intracellular, or secreted glycoproteins that are crucial in cellular processes, protein bioactivity and metabolic turnover. Modification by *O*-linked-*N*-acetyl glucosamine (*O*-GlcNAc) has rapidly emerged as a major cellular signaling mechanism with number of modified targets similar to protein phosphorylation. Many oncogenes and tumor suppressors are regulated by *O*-GlcNAcylation. *O*-GlcNAc modification is ubiquitous among eukaryotes, from yeast to humans and its modifying enzymes have been well characterized. *O*-GlcNAc modified nuclear and cytosolic targets include: transcription factors, signaling proteins, metabolic enzymes, mitochondrial trafficking, cell cycle regulation, glucose homeostasis. *O*-GlcNAc glycosylation is implicated in normal brain functions, etiology of neurodegeneration, type II diabetes, and pathways involved in morphogenesis and virulence factors of microbes and plant host cells. Since glycoproteins are not directly encoded in the genome, methods of characterization and analyses of glycoproteins are of great interest. Thus Assay Genie offers *O*-GlcNAc Modified Glycoprotein Assay Kit, a highly specific, simple and robust method for labeling and detection of *O*-GlcNAc-glycosylated proteins within cells. We use a modified glucosamine precursor that is fed directly into the cells, processed by the hexosamine pathway and incorporated into the proteins. Followed by click reaction with alkyne-containing dye, this system offers a powerful method for imaging the localization, trafficking, and dynamics of glycans, or detection by FACS for quantitative studies. Labeled Glycoproteins can be directly detected in 1D or 2D gels using the appropriate excitation sources, or enriched by immunoprecipitation with biotin-alkyne or antibodies prior to proteomic analysis. We provide sufficient materials for 100 assays in a 96-well plate format.

II. Applications:

- Identification and characterization of *O*-GlcNAc-glycosylated proteins within cells
- Imaging the localization, trafficking, and dynamics of modified glycans
- Detection and quantification of biosynthesis, subcellular localization and turnover of modified glycans
- Screening for genotoxic compounds and effectors of modified glycans in proliferating cells
- Evaluating effects of anti-cancer drugs and genotoxic agents on modified glycans

III. Sample Type:

- Suspension or adherent cell cultures

IV. Kit Contents:

Components	BN00934	Cap Code	Part Number
Wash Buffer (10X)	25 ml	NM	BN00934-1
Fixative Solution	10 ml	WM	BN00934-2
Permeabilization Buffer (10X)	25 ml	NM/Blue	BN00934-3
GlcAz Label (1000X)	10 µl	Clear	BN00934-4
Copper Reagent (100X)	100 µl	Blue	BN00934-5
Fluorescent Alkyne (100X)	100 µl	Green	BN00934-6
Reducing Agent (20X)	500 µl	Yellow	BN00934-7
Total DNA Stain (1000X)	20 µl	Blue	BN00934-8

V. User Supplied Reagents and Equipment:

- Tissue culture vessels and appropriate culturing media; flow cytometry vessels
- Phosphate Buffered Saline (PBS, pH 7.4)
- Sterile 0.1% Gelatin Solution (optional, only required for adhering suspension cells to the surface)
- Flow cytometer equipped with laser capable of excitation at 488 and 530/590 nm emission filters respectively
- Fluorescence microscope capable of excitation and emission at 440/490 and 540/580 nm respectively

VI. Storage Conditions and Reagent Preparation:

Upon arrival, store the entire kit at -20°C protected from light. Briefly centrifuge small vials prior to opening. Read the entire protocol before performing the assay.

- **10X Wash Buffer and 10X Permeabilization Buffers:** Dilute the 10X stocks 1:10 in sterile water, mix well. Store at 4 °C.
- **Fixative Solution:** Ready to use, after opening store at 4°C, protected from light.
- **Remaining components:** Store at -20°C protected from light. While in use, keep on ice and minimize light exposure.

VII. Assay Protocol:

Notes:

This assay was developed with HeLa (adherent) and Jurkat (suspension) cells and can be modified for any cell line. The protocol below refers to a 96-well tissue culture plate format and assay volume is 100 µl; adjust volumes accordingly for other plate formats. Growth conditions, cell number per well and other factors may affect the incorporation rate of the GlcAz Label; therefore optimize the assay for your cell type. We suggest an initial test of several GlcAz Label concentrations to find best conditions for your experimental design. Avoid stressing the cells by washes or temperature changes prior to incubation with GlcAz Label. All steps should be carried out at room temperature (RT) unless otherwise specified; equilibrate all buffers to RT prior to the experiment

1. Labeling with GlcAz Label:

- a. Seed the cell suspension of desired density directly into tissue culture vessels, or on coverslips for high resolution microscopy. **To immobilize suspension cells for microscopy:** add 100 μ l of 0.1% gelatin solution into each well of a tissue culture plate, tilt the plate to cover the entire well surface and place it in a tissue culture hood for 1 hour. Gently remove the gelatin solution and seed your cells. Allow the cells to recover overnight before the treatment.
- b. Next day, remove the media, and replace it with fresh aliquots containing 1X GlcAz Label. Include appropriate controls. **Negative control** -cells not exposed to the 1X GlcAz Label or treatment, **positive control** -cells incubated with 1X GlcAz Label only. *Do not add the GlcAz Label into the **negative control** cells.*
- c. Add treatments and incubate the cells for additional 1-3 days in a 37°C incubator, or for the period of time required by your experimental protocol. For analysis of trafficking and dynamics of cellular glycans take samples during incubation. Do not remove the drug-containing media while incubating with 1X GlcAz Label to avoid potential reversibility of drug action on label incorporation.
- d. Terminate the experiment, remove the media and rinse the cells once with 100 μ l of PBS, discard the supernatant. Always pellet the **suspension cells** at 300 x g for 5 min throughout the entire protocol. **For immobilized suspension cells:** Centrifuge the plate at 300 x g (or the lowest centrifuge setting) for 5 minutes to gently deposit the cells onto the surface. Tilt the plate and **gently** remove the media with a pipette tip. It is important to avoid excessive centrifugation speeds, which can damage the cells. *Make note of the place that is used, and perform subsequent aspirations from the same place.* Proceed to the Fixation and Permeabilization.

2. Fixation and Permeabilization:

For adherent and suspension cells: Add 100 μ l of Fixative Solution per well and incubate the cells for 15 min at RT protected from light. Remove the fixative and wash the cells twice with 100 μ l of 1X Wash Buffer. Remove the wash and add 100 μ l of 1X Permeabilization Buffer per well, incubate the cells for 10 min at RT. Remove the Permeabilization Buffer and replace it with a 20 μ l of fresh aliquot. Proceed to reaction and total DNA staining.

3. GlcAz reaction and total DNA staining:

- a. **Reaction cocktail:** Prepare 1X reaction cocktail according to the table below. Volumes should be multiplied by number of samples and reagents added in the exact order. Use the reaction cocktail within 15 minutes of preparation. *Cells should be protected from light during, and following the reaction and DNA staining.*

	Amount per Reaction
PBS	93 μ l
Copper Reagent (100X)	1 μ l
Fluorescent Alkyne (100X)	1 μ l
Reducing Agent (20X)	5 μ l

- b. **Reaction:** Add 100 μ l of 1X Reaction cocktail to each sample and incubate the cells for 30 min at room temperature protected from light. Remove the reaction cocktail and wash cells three times in 100 μ l of Wash Buffer. Remove the wash and suspend the cells in 100 μ l of PBS. Proceed to DNA staining. If no DNA staining is desired, proceed to Microscopic or FACS analysis. **DNA staining:** Prepare 1X dilution of Total DNA Stain and add 100 μ l per well. Incubate the cells for 20 minutes at room temperature, or refrigerate at 4 °C protected from light. Remove the stain solution and replace with 100 μ l of PBS.

4. **Fluorescence Microscope analysis:** Examine labeled glycoproteins using FITC filter and UV laser for total DNA staining. **FACS analysis:** Transfer the 100 μ l cell suspension into flow cytometry vessels. Analyze samples in FL-1 channel for signal generated by labeled glycoproteins.




Figure: Analysis of metabolic labeling of GlcAz labeled glycans in proliferating cells. **(A)** Jurkat cells (1×10^6 cells/ml) were cultured in presence of 1X GlcAz Label for 24 hours at 37°C. Modified glycoproteins were detected according to the kit protocol and green fluorescence was analyzed by FACS (FL-1 channel). Negative control (**white line**), Background control (**purple line**), fluorescence corresponding to intracellular O-GlcNAc-glycosylated proteins (**green line**). **(B)** Fluorescence Microscope analysis of intracellular O-GlcNAc-glycosylated proteins in HeLa cells. High resolution images (middle and right panels) clearly show cytoplasmic and nuclear localization of GlcAz modified glycans.

FOR RESEARCH USE ONLY! Not to be used on humans.