

Product Datasheet **GenieFluor 647 Anti-Human CD49d Antibody [BU49]** Catalogue Code: AGEL0888

Antibody Data

Product SKU:	AGEL0888	Clone:	BU49
Applications:	FCM		
Reactivity:	Human		

Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

Product Information:

Alternate Names: Uniprot ID: Background:	Integrin alpha-4;ITGA4;CD49 antigen-like family member D;Integrin alpha-IV;VLA-4 subunit alpha;CD49d; P13612 CD49d is a 150 kD α integrin chain known as α 4 integrin or VLA-4 α chain. It forms a heterodimer with either integrin β 1 (α 4 β 1, VLA-4) or β 7 (α 4 β 7). CD49d is expressed broadly on T lymphocytes, B lymphocytes, monocytes, thymocytes, eosinophils, basophils, mast cells, NK cells, dendritic cells, and some non-hematopoietic cells, but not on normal red blood cells, platelets or neutrophils. VLA-4 binds to VCAM-1 (CD106) and fibronectin. α 4 β 7 is the receptor for VCAM-1 and MAdCAM-1. CD49d participates in mononuclear cell trafficking to endothelial sites of inflammation and has roles in cell-cell interactions and cell adhesion to extracellular matrices. CD49d is involved in lymphocyte migration, T cell activation, and hematopoietic stem cell differentiation. CD49d is a marker to isolate pure populations of Treg cells due to its absence on Foxp3+ cells.		
Form:	Liquid	647 Excitation and Emission Spectra	
Conjugation:	Genie Fluor647	100 -	
Size:	20 Tests, 100 Tests, 200 Tests	80 · (2)	
Host Species:	Mouse	W 40 .	
Isotype:	Mouse IgG1, к	E 40 0 0 0 0 0 0 0 0 0 0 0 0 0	
Isotype Control:	Genie Fluor 647 Mouse IgG1, κ Isotype Control[MOPC-21] [Product AGEL0888]		
Storage Buffer:	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.		
Shipping:	Biological ice pack at 4°C		



- **Stability & Storage:** Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.
- **Recommended** Usage: Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μ L staining volume or per 100 μ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.