

Product Datasheet **PE/Cyanine7 Anti-Human CD117 Antibody [104D2]** Catalogue Code: AGEL3067

Antibody Data

Product SKU:	AGEL3067	Clone:	104D2
Applications:	FCM		
Reactivity:	Human		

Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

Product Information:

Alternate Names: Uniprot ID: Background:	Mast/stem cell growth factor receptor Kit;Kit;SCFR;Proto-oncogene c-Kit;Tyrosine-protein kinase Kit;CD117; P10721 CD117 is a 145 kD protein tyrosine kinase also known as c-Kit. It is a receptor for stem cell factor or c-Kit ligand. CD117 is expressed on pluripotent hematopoietic progenitor cells (approximately 1-4% bone marrow cells), mast cells, and acute myeloid leukemia cells (AML). CD117 binding of c-Kit ligand induces phosphorylation of CD117 and stimulates proliferation and survival of primitive hematopoietic stem cells as well as erythroid- committed and granulo-monocytic committed cells.		
Form:	Liquid	PE/Cyanine7 Excitation and Emission Spectra	
Conjugation:	PE/Cyanine 7	100 -	
Size:	20 Tests, 100 Tests, 200 Tests		
Host Species:	Mouse	\$ 60 - 1 - 1 - 1 - 1 - 1 - 1 - 1 - 1 - 1 -	
Isotype:	Mouse IgG1, κ	20 0 350 400 450 500 550 600 650 700 750 800 850 Wavelength (nm)	
		Ex:495;565;755 nm; Em:775 nm	

Isotype Control: PE/Cyanine7 Mouse IgG1, κ Isotype Control[MOPC-21] [Product AGEL3067]

Storage Buffer: Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.

Shipping: Biological ice pack at 4°C



- **Stability & Storage:** Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.
- **Recommended** Usage: Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μ L staining volume or per 100 μ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.