

Product Datasheet GenieFluor Violet 450 Anti-Human CD20 Antibody [2H7] Catalogue Code: AGEL3090

Antibody Data

Product SKU:	AGEL3090	Clone:	2H7
Applications:	FCM		
Reactivity:	Human		

Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

Product Information:

Alternate Names:	Bp35;Leukocyte surface antigen Leu-16;MS4A1;B-lymphocyte surface antigen B1;		
Uniprot ID:	P11836		
Background:	CD20 is a 33-37 kD, four transmembrane spanning protein, also known as B1 and Bp35. CD20 is expressed on pre-B-cells, resting and activated B cells (not plasma cells), some follicular dendritic cells, and at low levels on a T cell subset. CD20 is heavily phosphorylated on activated B cells and malignant B cells. Homo-oligomeric complexes of CD20 are thought to form Ca2+ conductive ion channels in the plasma membrane of B cells. The CD20 molecule is involved in B-cell activation and is associated with various Src family kinases (Lyn, Lck, Fyn). It exists in a complex with MHC class I and II, CD53, CD81, and CD82.		
Form:	Liquid	450 Excitation and Emission Spectra	
Conjugation:	Genie FluorViolet 450	100	
Size:	20 Tests, 100 Tests, 200 Tests	80	
Host Species:	Mouse	(%) 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.	
Isotype:	Mouse IgG2b, к	$\frac{1}{20} \frac{1}{\frac{1}{350} \frac{1}{400} \frac{1}{450} \frac{1}{500} \frac{1}{550} \frac{1}{550} \frac{1}{600} \frac{1}{650} \frac{1}{700}}{\frac{1}{100} \frac{1}{100} \frac{1}{100}$	
Isotype Control:	Genie Fluor Violet 450 Mouse IgG2b, κ Isotype Control[MPC-11] [Product AGEL3090]		

Storage Buffer: Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.

Shipping: Biological ice pack at 4°C



- **Stability & Storage:** Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.
- **Recommended** Usage: Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μ L staining volume or per 100 μ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.