

Product Datasheet **PE/GenieFluor 594 Anti-Human CD45RO Antibody [UCHL1]** Catalogue Code: AGEL3148

## Antibody Data

Product SKU:	AGEL3148	Clone:	UCHL1
Applications:	FCM		
Reactivity:	Human		

## Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

## **Product Information:**

Alternate Names:	-		
Uniprot ID:	-		
Background:	CD45RO is a 180 kD single chain membrane glycoprotein. It is a splice variant of tyrosine phosphatase CD45, lacking the A, B, and C determinants. The CD45RO isoform is expressed on activated and memory T cells, some B cell subsets, activated monocytes/macrophages, and granulocytes. CD45RO enhances both T cell receptor and B cell receptor signaling mediated activation. CD45 and its isoforms non-covalently associate with lymphocyte phosphatase-associated phosphoprotein (LPAP) on T and B lymphocytes. CD45 has been reported to be associated with several other cell surface antigens including CD1, CD2, CD3, and CD4. CD45 has also been reported to bind galectin-1 and CD22. CD45 isoform expression can change in response to cytokines.		
Form:	Liquid	594 Excitation and Emission Spectra	
Conjugation:	PE/Genie Fluor594	100	
Size:	20 Tests, 100 Tests, 200 Tests	80 - Se	
Host Species:	Mouse	(%) pe 60	
Isotype:	Mouse IgG2a, к	E 40 20 0 350 400 450 500 550 600 650 700 750 800 850 Wavelength (nm)	
Isotype Control:	PE/Genie Fluor 594 Mouse IgG2a, κ Isotype Control[C1.18.4] [Product AGEL3148]		
Storage Buffer:	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.		

Shipping: Biological ice pack at 4°C



- **Stability & Storage:** Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.
- **Recommended** Usage: Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5  $\mu$ L of antibody per test (million cells in 100  $\mu$ L staining volume or per 100  $\mu$ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.