



フローサイトメトリー用抗体



13のカラフルな
蛍光色素！

3,500超のフローサイトメトリー検証済み抗体からお選びください！

抗ヒト | 抗マウス | 抗ラット | 抗非ヒト霊長類

13の蛍光色素 | 20~200テストパックと25 ug~1 mgパック

ミッションとブランドプロミス

ミッション

私たちのミッションは、革新的な製品、当社チームと小さな
精霊との優れた関係を通じて、科学者を笑顔にすることです。



COLM RYAN PhD
CEO & co-founder of
Assay Genie



SEÁN MAC FHEARRAIGH
PhD CTO & co-founder of
Assay Genie

最高のサポートと保証

Assay Genieは全てのお客様に最高のサポートを提供します！

実験前後のアプリケーションベースのサポートだけでなく、稀に問題が生じた場合でも、決まった「顧客中心」のステップで、お客様にご満足いただけるよう努めています。

ですので、心配なさらないでください。当社の製品が指定通りにワークしない場合でも、100%の返金保証を提供しています。



テクニカルサポートメール、動画、電話

トラブルシューティング

交換品のキット

100%返金保証

世界中へ迅速な納品

信頼のローカルディストリビューターから購入する場合でも、当社の直接購入ネットワークの場合でも、常に時間通りに製品をお届けするように努めています。


www.info@assaygenie.comから24時間年中無休でお問い合わせいただき、お客様の研究室へのお届け時間をご確認ください。



フローサイトメトリー用抗体

Assay Genieの優れたフローサイトメトリー用抗体で、研究の可能性を解き放ちます。効率と精度を重視した設計の当社製品は、信頼性と品質の高い研究に理想的な選択肢です。
当社抗体独自の特長をご覧ください。

- 多様な蛍光色素の選択肢：**13の異なる蛍光色素を提供。当社製品は様々なパネルデザインを可能にします。
- 幅広いターゲットの範囲：**35,000を超える細胞内・細胞外タンパク質のターゲットからお選びください。
- 検証済みの標識と親和性：**各抗体は最適な標識と強い親和性を保証するため、厳格に試験されています。
- 最上級のモノクローナル抗体：**当社のフローサイトメトリー検証済み抗体は最高品質。お客様の実験で正確性と一貫性を保証しています。
- 最先端の標識技術：**先端技術を利用し、当社の抗体は免疫細胞の効率的な同定、および細胞活性化状態の正確な評価を容易にします。
- 研究に合わせカスタマイズ：**当社のコレクションからお探しいただき、お客様のフローサイトメトリー研究要件に最適な製品を見つけてください。



フローサイトメトリー用蛍光色素

蛍光色素	明るさ	蛍光色	レーザー	最大励起波長（nm）	最大蛍光波長（nm）
GenieFluor Violet 450	✓ ✓	Blue	405	410	450
PE	✓ ✓ ✓ ✓ ✓	Yellow	488, 561	495, 565	575
FITC	✓	Green	488	490	530
GenieFluor 488	✓ ✓ ✓	Green	488	595	520
PerCP	✓	Red	488	440, 480, 675	675
PE/Cyanine5	✓ ✓ ✓ ✓ ✓	Red	488, 561	495, 565, 655	670
PE/Cyanine5.5	✓ ✓ ✓ ✓	Red	488, 561	495, 565, 675	690
PE/GenieFluor 594	✓ ✓ ✓ ✓ ✓	Red	488, 561	566	610
APC	✓ ✓ ✓ ✓	Red	633	650	660
GenieFluor 647	✓ ✓ ✓ ✓	Red	633	650	670
PerCP/Cyanine5.5	✓ ✓ ✓	Far Red	488	440, 480, 675	675
GenieFluor Red 780	✓ ✓	Infrared	633	625	765
PE/Cyanine7	✓ ✓ ✓ ✓ ✓	Infrared	488, 561	495, 565, 755	775

フローサイトメトリー用抗体 | 抗ヒト

名称	反応性	クローン	GeniePure	未標識	Biotin	FITC	PE	APC	PerCP	PE/Cy5	PE/Cy7	PE/Cy5.5	PerCP/Cy5.5	GenieFluor 488	GenieFluor 647	PE/GenieFluor 594	GenieFluor Red 780	GenieFluor Violet 450
CD1a	Anti-Human	OKT-6	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		
CD2	Anti-Human	TS1/8				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
CD3	Anti-Human	UCHT1	✓	✓	✓	✓	✓	✓	✓		✓		✓	✓			✓	✓
CD3e	Anti-Human	OKT-3	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD4	Anti-Human	RPA-T4	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD5	Anti-Human	UCHT2		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
CD7	Anti-Human	124-1D1		✓		✓	✓	✓								✓		✓
CD8a	Anti-Human	OKT-8		✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓		✓	✓
CD8a	Anti-Human	HIT8a		✓		✓	✓							✓				
CD9	Anti-Human	HI9a		✓	✓	✓	✓	✓		✓	✓	✓						✓
CD10	Anti-Human	CB-CALLA		✓	✓	✓	✓	✓		✓	✓			✓	✓			
CD10	Anti-Human	HI10a		✓		✓	✓	✓		✓	✓			✓	✓			
CD11a	Anti-Human	TS1/22.1.1.13		✓	✓	✓	✓	✓						✓	✓			
CD11a	Anti-Human	R7-1	✓	✓		✓	✓	✓		✓				✓	✓	✓		
CD11b	Anti-Mouse/Human	M1/70	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD11b	Anti-Human	ICRF44				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	
CD11c	Anti-Human	BU15				✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		
CD13	Anti-Human	B-F10				✓	✓	✓			✓		✓	✓		✓		
CD14	Anti-Human	M5E2				✓	✓	✓		✓	✓	✓	✓	✓		✓	✓	✓
CD15	Anti-Human	HI98		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓				✓
CD15	Anti-Human	W6D3				✓	✓	✓		✓	✓	✓	✓	✓	✓			
CD16	Anti-Human	3G8		✓		✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	
CD18	Anti-Human	TS1/18.1.2.11		✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		✓
CD19	Anti-Human	CB19			✓	✓	✓	✓			✓			✓		✓		✓
CD19	Anti-Human	HI19a		✓		✓	✓	✓		✓	✓				✓	✓		✓
CD20	Anti-Human	BCA/B20				✓	✓	✓	✓				✓	✓	✓	✓	✓	✓
CD20	Anti-Human	2H7	✓	✓		✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
CD20/PD-1	Anti-Human	J110	✓	✓														
CD21	Anti-Human	BU32		✓	✓	✓	✓	✓		✓	✓		✓	✓	✓			
CD22	Anti-Human	BU59				✓	✓	✓		✓			✓	✓		✓		✓
CD22	Anti-Human	HIB22		✓		✓	✓	✓										
CD24	Anti-Human	ML5				✓	✓	✓		✓	✓	✓	✓	✓				
CD25	Anti-Human	BC96				✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		
CD27	Anti-Human	O323				✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	
CD28	Anti-Human	CD28.2	✓	✓		✓	✓	✓		✓	✓		✓	✓	✓	✓		
CD29	Anti-Human	TS2/16.2.1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
CD31	Anti-Human	158-2B3				✓	✓	✓		✓			✓	✓	✓			
CD32	Anti-Human	IV-3		✓	✓	✓	✓	✓		✓	✓	✓		✓	✓	✓		✓
CD326	Anti-Human	9C4				✓	✓	✓		✓	✓							✓
CD33	Anti-Human	6C5				✓	✓	✓			✓	✓		✓				✓
CD34	Anti-Human	581				✓	✓	✓					✓			✓		
CD34	Anti-Human	4H11		✓		✓				✓				✓				
CD35	Anti-Human	E11		✓	✓	✓	✓			✓	✓			✓	✓	✓		
CD36	Anti-Human	5-271				✓	✓							✓				

フローサイトメトリー用抗体

名称	反応性	クローン	GeniePure	未標識	Biotin	FITC	PE	APC	PerCP	PE/Cy5	PE/Cy7	PE/Cy5.5	PerCP/Cy5.5	GenieFluor 488	GenieFluor 647	PE/GenieFluor 594	GenieFluor Red 780	GenieFluor Violet 450
CD37	Anti-Human	IPO-24				✓	✓						✓	✓	✓	✓		
CD38	Anti-Human	HIT2		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
CD39	Anti-Human	A1				✓	✓	✓		✓	✓	✓		✓	✓		✓	
CD40	Anti-Human	G28.5	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓
CD41	Anti-Human	HIP8		✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
CD44	Anti-Human	P2A1				✓	✓	✓			✓		✓	✓	✓			
CD44	Anti-Human/Mouse	IM7	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓
CD44	Anti-Human	Heremes-1	✓	✓		✓	✓	✓		✓	✓		✓	✓	✓			✓
CD45	Anti-Human	HI30		✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD45RA	Anti-Human	HI100		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓
CD45RO	Anti-Human	UCHL1		✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
CD47	Anti-Human/Ms/Rat	MIAP410	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓				✓
CD47	Anti-Human	CC2C6D4		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
CD48	Anti-Human	156-4H9				✓	✓	✓		✓		✓				✓		✓
CD49d	Anti-Human	BU49				✓									✓			
CD49d	Anti-Human	9F10				✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		
CD49d	Anti-Human	HP1/2				✓	✓				✓		✓					
CD54	Anti-Human	15.2				✓	✓	✓			✓		✓			✓		
CD56	Anti-Human	5.1H11				✓	✓	✓									✓	✓
CD56	Anti-Human	MY31				✓												
CD56	Anti-Human	B-A19		✓		✓	✓	✓			✓	✓			✓			
CD57	Anti-Human	HNK-1		✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
CD58	Anti-Human	TS2/9/1		✓	✓		✓	✓		✓	✓			✓		✓		
CD61	Anti-Human	VI-PL2				✓	✓	✓				✓				✓		
CD62L	Anti-Human	DREG56		✓	✓	✓	✓	✓		✓	✓		✓	✓	✓	✓	✓	✓
CD62L	Anti-Human	HI62L				✓	✓	✓		✓								
CD64	Anti-Human	10.1		✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓
CD66b	Anti-Human	G10F5				✓	✓	✓			✓						✓	✓
CD68	Anti-Human	Y1/82A																
CD69	Anti-Human	FN50				✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
CD73	Anti-Human	AD2				✓	✓	✓		✓	✓	✓		✓	✓		✓	
CD74	Anti-Human	LN2		✓		✓	✓	✓								✓		✓
CD80	Anti-Human	2D10					✓	✓			✓							
CD81	Anti-Human	1.3.3.22				✓	✓	✓						✓		✓		
CD83	Anti-Human	HB15e					✓	✓			✓			✓		✓		
CD86	Anti-Human	BU63		✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	
CD90	Anti-Human	5E10				✓	✓	✓		✓	✓	✓	✓		✓	✓		
CD95	Anti-Human	DX2				✓	✓	✓		✓	✓			✓	✓			
CD105	Anti-Human	SN6					✓	✓							✓			✓
CD107a/LAMP-1	Anti-Human	H4A3				✓	✓	✓					✓	✓		✓		
CD117	Anti-Human	104D2				✓	✓	✓			✓							
CD122/IL-2RB	Anti-Human	TU27				✓	✓	✓										
CD123	Anti-Human	6H6			✓		✓			✓			✓	✓	✓			
CD127/IL-7RA	Anti-Human	A019D5			✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓
CD133	Anti-Human	W6B3C1					✓	✓										

フローサイトメトリー用抗体 | 抗ヒト・抗ラット

名称	反応性	クローン	GeniePure	未標識	Biotin	FITC	PE	APC	PerCP	PE/Cy5	PE/Cy7	PE/Cy5.5	PerCP/Cy5.5	GenieFluor 488	GenieFluor 647	PE/GenieFluor 594	GenieFluor Red 780	GenieFluor Violet 450
CD134	Anti-Human	Ber-ACT35				✓	✓	✓									✓	
CD138	Anti-Human	DL-101						✓										
CD147	Anti-Human	HIM6		✓	✓	✓	✓	✓		✓	✓		✓	✓	✓	✓		
CD154/CD40L	Anti-Human/Monkey	5C8	✓	✓	✓													
CD161	Anti-Human	HP-3G10				✓	✓	✓		✓	✓	✓	✓	✓	✓			
CD183/CXCR3	Anti-Human	G025H7				✓	✓	✓			✓							
CD184/CXC4	Anti-Human	12G5					✓	✓										
CD196/CCR6	Anti-Human	CCR6				✓	✓	✓			✓							
CD196/CCR6	Anti-Human	G034E3				✓	✓	✓			✓				✓			✓
CD197/CCR7	Anti-Human	G043H7				✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓
CD200	Anti-Human	OX-104					✓	✓			✓							
CD206	Anti-Human	15_2						✓										
CD235	Anti-Human	HIR2		✓	✓	✓	✓	✓		✓	✓			✓	✓			
CD273/PD-L2	Anti-Human	24F.10C12					✓	✓										
CD279/PD-1	Anti-Human	J116	✓	✓														
CD279/PD-1	Anti-Human	EH12.2H7				✓	✓	✓					✓					
CD303	Anti-Human	201A				✓		✓										
HLA-A,B,C	Anti-Human	W6/32	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			✓
HLA-DR	Anti-Human	L243	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓				✓	✓
IFN-γ	Anti-Human	B27	✓	✓		✓	✓	✓										
IgD	Anti-Human	IA6-2				✓	✓	✓						✓				✓
IgM	Anti-Human	MHM-88				✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
IL-10	Anti-Human	JES3-9D7				✓	✓	✓						✓	✓			
IL-17A	Anti-Human	BL168				✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
IL-2	Anti-Human	MQ1-17H12				✓	✓	✓										
IL-21	Anti-Human	3A3-N2					✓	✓										
IL-4	Anti-Human	MP4-25D2	✓	✓		✓	✓	✓								✓		
IL-5	Human/Mouse	TRFK5	✓	✓			✓	✓										
IL-6	Anti-Human	MQ2-13A5				✓	✓	✓										
KLRG-1	Anti-Human/Mouse	2F1		✓		✓	✓	✓			✓			✓				
TCR γ/δ	Anti-Human	B1				✓												

名称	反応性	クローン	GeniePure	未標識	Biotin	FITC	PE	APC	PerCP	PE/Cy5	PE/Cy7	PE/Cy5.5	PerCP/Cy5.5	GenieFluor 488	GenieFluor 647	PE/GenieFluor 594	GenieFluor Red 780	GenieFluor Violet 450
CD3	Anti-Rat	G4.18				✓	✓	✓					✓			✓		
CD4 (domain 1)	Anti-Rat	OX-38		✓	✓	✓	✓	✓	✓	✓	✓		✓	✓			✓	✓
CD44H	Anti-Rat	OX-49				✓	✓	✓		✓	✓		✓					
CD45	Anti-Rat	OX-1				✓	✓	✓		✓	✓		✓				✓	✓
CD47	Anti-Human/Ms/Rat	MIAP410	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓				✓
CD8a	Anti-Rat	OX-8				✓	✓	✓			✓	✓		✓				
CD90	Anti-Rat / Mouse	OX-7			✓	✓	✓	✓			✓		✓				✓	

フローサイトメトリー用抗体 | 抗マウス

名称	反応性	クローン	GeniePure	未標識	Biotin	FITC	PE	APC	PerCP	PE/Cy5	PE/Cy7	PE/Cy5.5	PerCP/Cy5.5	GenieFluor 488	GenieFluor 647	PE/GenieFluor 594	GenieFluor Red 780	GenieFluor Violet 450
CD3	Anti-Mouse	17A2	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD1d	Anti-Mouse	19G11	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			
CD3e	Anti-Mouse	145-2C11	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
CD4	Anti-Mouse	GK1.5	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD4	Anti-Mouse	RM4-5					✓	✓			✓		✓				✓	✓
CD5	Anti-Mouse	53-7.3			✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		✓	✓
CD11a	Anti-Mouse	FD441.8	✓	✓	✓	✓	✓	✓		✓	✓	✓		✓	✓	✓		
CD11b	Anti-Mouse/Human	M1/70	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD11c	Anti-Mouse	N418			✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓
CD14	Anti-Mouse	Sa14-2				✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓
CD16/32	Anti-Mouse	2.4G2	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓
CD18	Anti-Mouse	M18/2	✓	✓	✓	✓	✓				✓			✓				
CD19	Anti-Mouse	1D3	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓
CD22	Anti-Mouse	Cy34.1	✓	✓	✓			✓		✓	✓	✓	✓	✓	✓			✓
CD23	Anti-Mouse	B3B4				✓	✓	✓			✓			✓	✓			
CD24	Anti-Mouse	M1/69				✓	✓	✓		✓	✓	✓	✓	✓	✓			
CD25	Anti-Mouse	PC-61.5.3	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓
CD28	Anti-Mouse	37.51	✓	✓			✓	✓						✓	✓			
CD31	Anti-Mouse	390			✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			
CD34	Anti-Mouse	RAM34																
CD36	Anti-Mouse	HM36																
CD38	Anti-Mouse	NIM45	✓	✓		✓	✓	✓			✓							
CD40	Anti-Mouse	FGK4.5/FGK45	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		
CD41	Anti-Mouse	MWReg30				✓	✓	✓			✓							
CD44	Anti-Human/Mouse	IM7	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓
CD45	Anti-Mouse	30-F11				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD45.1	Anti-Mouse	A20				✓	✓	✓										✓
CD45.2	Anti-Mouse	104.2	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
CD45R/B220	Anti-Mouse	RA3.3A 1/6.1	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD47	Anti-Human/Ms/Rat	MIAP410	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓				✓
CD48	Anti-Mouse	HM48-1				✓	✓	✓	✓		✓		✓	✓	✓		✓	
CD49b	Anti-Mouse	DX5			✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓
CD51	Anti-Mouse	RMV-7			✓		✓	✓			✓							
CD54	Anti-Mouse	YN/1.7.4	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			
CD62L	Anti-Mouse	Mel14	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			
CD64	Anti-Mouse	X54-5F/7.1				✓	✓	✓		✓	✓						✓	
CD69	Anti-Mouse	H1.2F3				✓	✓	✓			✓						✓	
CD71	Anti-Mouse	R17 217.1.3/TIB-219	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		✓
CD73	Anti-Mouse	TY/23	✓	✓	✓	✓	✓	✓			✓				✓			
CD80	Anti-Mouse	16-10A1	✓	✓	✓	✓	✓	✓		✓	✓	✓				✓	✓	✓
CD86	Anti-Mouse	GL-1	✓	✓	✓	✓	✓	✓			✓	✓		✓	✓	✓		
CD8a	Anti-Mouse	53-6.7	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CD90	Anti-Rat / Mouse	OX-7			✓	✓	✓	✓			✓		✓				✓	
CD90	Anti-Mouse	M5/49.4.1	✓	✓	✓	✓	✓	✓			✓							
CD90.1	Anti-Mouse	1.90E+13	✓	✓														
CD90.2	Anti-Mouse	30H12	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			✓
CD103	Anti-Mouse	M290	✓	✓	✓	✓	✓	✓		✓	✓	✓		✓	✓			
CD105	Anti-Mouse	MJ7/18				✓	✓	✓										
CD106	Anti-Mouse	M/K-2.7	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		

フローサイトメトリー用抗体 | 抗マウス

名称	反応性	クローン	GeniePure	未標識	Biotin	FITC	PE	APC	PerCP	PE/Cy5	PE/Cy7	PE/Cy5.5	PerCP/Cy5.5	GenieFluor 488	GenieFluor 647	PE/GenieFluor 594	GenieFluor Red 780	GenieFluor Violet 450
CD107a	Anti-Mouse	1D4B				✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		✓
CD115	Anti-Mouse	AFS98	✓	✓			✓	✓										
CD117	Anti-Mouse	2B8	✓	✓	✓	✓	✓	✓			✓			✓	✓			
CD119	Anti-Mouse	GR-20	✓	✓	✓		✓	✓										
CD120b	Anti-Mouse	TR75-54.7	✓	✓	✓													
CD122	Anti-Mouse	5H4	✓	✓	✓	✓	✓	✓		✓	✓	✓		✓	✓			
CD127	Anti-Mouse	A7R34		✓	✓	✓	✓	✓		✓	✓		✓	✓	✓			
CD132	Anti-Mouse	3 E12	✓	✓	✓													
CD150	Anti-Mouse	TC15-12F12.2				✓	✓	✓		✓	✓	✓		✓	✓			
CD162	Anti-Mouse	4RA10	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		✓
CD183/CXCR3	Anti-Mouse	CXCR3-173	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			
CD200	Anti-Mouse	OX-90				✓		✓			✓							
CD206	Anti-Mouse	C068C2				✓	✓	✓			✓							
CD209b	Anti-Mouse	22D1	✓	✓	✓													
CD210	Anti-Mouse	1B1.3A	✓	✓	✓													
CD272	Anti-Mouse	PK18.6	✓	✓	✓	✓	✓	✓			✓							
CD274/PD-L1	Anti-Mouse	10F.9G2	✓	✓	✓	✓	✓	✓				✓		✓	✓			
CD274/PD-L1	Anti-Mouse	29E.2A3	✓	✓	✓	✓	✓	✓		✓	✓		✓					
CD279/PD-1	Anti-Mouse	29F.1A12	✓	✓	✓	✓	✓	✓			✓	✓	✓		✓		✓	
CD326	Anti-Mouse	G8.8				✓	✓	✓			✓			✓				
CD335	Anti-Mouse	29A1.4				✓	✓	✓		✓	✓	✓				✓		
CD366/Tim-3	Anti-Mouse	RMTM-23	✓	✓		✓	✓	✓										
F4/80	Anti-Mouse	Cl:A3-1	✓	✓	✓	✓	✓	✓		✓	✓	✓		✓	✓	✓	✓	✓
FcεRIα	Anti-Mouse	MAR.1				✓	✓	✓										
Foxp3	Anti-Mouse	3G3				✓	✓	✓										
H-2	Anti-Mouse	M1/42	✓	✓	✓	✓	✓	✓		✓		✓	✓	✓	✓		✓	
IFN-γ	Anti-Mouse	XMG1.2	✓	✓	✓	✓	✓	✓								✓		
IgD	Anti-Mouse	11-26c.2a			✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	
IgM	Anti-Mouse	RMM-1			✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			
IL-10	Anti-Mouse	JES5-16E3				✓	✓	✓										
IL-17A	Anti-Mouse	TC11-18H10.1			✓	✓	✓	✓							✓			
IL-17A	Anti-Mouse	17F3		✓		✓	✓	✓			✓							✓
IL-2	Anti-Mouse	JES6-5H4	✓	✓		✓	✓	✓										
IL-4	Anti-Mouse	11B11	✓	✓			✓	✓										
IL-5	Human/Mouse	TRFK5	✓	✓			✓	✓										
IL-6	Anti-Mouse	MP5-20F3	✓	✓			✓	✓										
KLRG-1	Anti-Human/Mouse	2F1		✓		✓	✓	✓			✓			✓				
Ly6A/E(Sca-1)	Anti-Mouse	D7				✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
Ly6C	Anti-Mouse	Monts 1	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓
Ly6G	Anti-Mouse	1A8	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			✓
Ly-6G/Ly-6C (Gr-1)	Anti-Mouse	RB6-8C5	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓
MHC-II	Anti-Mouse	M5/114	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			✓
MPO	Anti-Mouse	ABT-101007				✓	✓	✓							✓			
NK1.1	Anti-Mouse	PK136	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			
Perforin	Anti-Mouse	S16009A																
TCR γ/δ	Anti-Mouse	UC7-13D5	✓	✓														
TCR γ/δ	Anti-Mouse	GL3			✓	✓												
TCRβ	Anti-Mouse	H57-597 (HB218)	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓
TER-119	Anti-Mouse	TER-119	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓		

アイソタイプコントロールとカクテル

名称	反応性	クローン	GeniePure	未標識	Biotin	FITC	PE	APC	PerCP	PE/Cy5	PE/Cy7	PE/Cy5.5	PerCP/Cy5.5	GenieFluor 488	GenieFluor 647	PE/GenieFluor 594	GenieFluor Red 780	GenieFluor Violet 450
Armenian Hamster IgG Isotype Control	Isotype Control	PIP	✓	✓	✓	✓	✓	✓			✓	✓	✓	✓	✓	✓		✓
Mouse IgG1, κ Isotype Control	Isotype Control	MOPC-21	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓
Mouse IgG2a, κ Isotype Control	Isotype Control	C1.18.4	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓
Mouse IgG2b, κ Isotype Control	Isotype Control	MPC-11	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓
Mouse IgG3, κ Isotype Control	Isotype Control	A112-3		✓	✓	✓	✓	✓		✓			✓	✓	✓	✓	✓	✓
Mouse IgM, κ Isotype Control	Isotype Control	MM-30		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rat IgG1, κ Isotype Control	Isotype Control	HRPN	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓
Rat IgG2a, κ Isotype Control	Isotype Control	2A3	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rat IgG2a, κ Isotype Control	Isotype Control	RTK2758														✓		
Rat IgG2a, λ Isotype Control	Isotype Control	B-39-4	✓		✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓
Rat IgG2a, λ Isotype Control	Isotype Control	G013C12														✓		
Rat IgG2b, κ Isotype Control	Isotype Control	LTF-2	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rat IgM, κ Isotype Control	Isotype Control	RTK2118	✓		✓	✓	✓	✓		✓				✓	✓	✓		✓
Syrian Hamster IgG Isotype Control	Isotype Control	SHG-1	✓		✓	✓	✓							✓	✓			✓

名称	反応性	クローン	FITC	PE	APC	PerCP	PE/Cy5	PE/Cy7	PE/Cy5.5	PerCP/Cy5.5
CD3-FITC/CD19-APC/CD16/CD56-PE Cocktail	Anti-Human	FITC CD3 [OKT-3], APC CD19 [4G7], PE CD16 [3G8], PE CD56 [5.1H11]	✓	✓	✓					
CD19-FITC/CD56-PE/CD3-PE/Cy7/CD45-PerCP Cocktail	Anti-Human	FITC CD19 [CB19], PE CD56 [5.1H11], PE/Cyanine7 CD3 Antibody [OKT-3], PerCP CD45 [HI30]	✓	✓		✓		✓		
CD3-APC/CD4-FITC/CD8a-PE Cocktail	Anti-Human	APC CD3 [OKT-3], FITC CD4 [RPA-T4], PE CD8a [OKT-8]	✓	✓	✓					
CD3-APC/CD4-FITC/CD8a-PerCP-Cy5.5 Cocktail	Anti-Human	APC CD3 [OKT-3], FITC CD4 [RPA-T4], PerCP/Cyanine5.5 CD8a [OKT-8]	✓		✓				✓	✓
CD3-FITC/CD16+CD56-PE Cocktail	Anti-Human	FITC CD3 [OKT-3], PE CD16 [3G8], PE CD56 [5.1H11]	✓	✓						
CD3-FITC/CD19-PE Cocktail	Anti-Human	FITC CD3 [OKT-3], PE CD19 Antibody [CB19]	✓	✓						
CD3-FITC/CD4-PE/CD45-PerCP-Cy5.5 Cocktail	Anti-Human	FITC CD3 [OKT-3],PE CD4 [RPA-T4], PerCP/Cyanine5.5 CD45 [HI30]	✓	✓						✓
CD3-FITC/CD4-PE/CD8a-PE/Cy/CD45-PerCP Cocktail	Anti-Human	FITC CD3 [OKT-3], PE CD4 [RPA-T4], PE/Cyanine7 CD8a [OKT-8], PerCP CD45 [HI30]	✓	✓		✓		✓		
CD8a-FITC/CD4-PE/CD3-PE/Cy7/CD45-PerCP Cocktail	Anti-Human	FITC CD8a Antibody [OKT-8], PE CD4 Antibody [RPA-T4], PE/Cyanine7 CD3 [OKT-3], PerCP CD45 [HI30]	✓	✓		✓		✓		



フィードバックをお待ちしています！

当社の成長とアイデアは、笑顔の科学者に支えられています。
当社のコミュニティと会話にご参加ください！

フローサイトメトリー用抗体は、www.assaygenie.com
または当社の世界的なディストリビューターネットワーク
よりご注文ください。

www.assaygenie.com | Assay Genie 2022

25 Windsor Place, Dublin, Ireland | 2 Duke St., London, UK

