



Recombinant Protein Technical Manual

Recombinant Rat IL13RA2/CD213A2 Protein (Fc Tag)(Active)
RPES0504

Product Data:

Product SKU: RPES0504

Size: 100µg

Species: Rat

Expression host: HEK293 Cells

Uniprot: Q8VHK6

Protein Information:

Molecular Mass: 63.4 kDa

AP Molecular Mass: 89-93 kDa

Tag: C-Fc

Bio-activity: Measured by its ability to inhibit IL13-dependent proliferation of TF human erythroleukemic cells. The ED50 for this effect is typically 0.005-0.03 µg/mL.

Purity: > 90 % as determined by SDS-PAGE

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: IL13R;IL13RA2

Immunogen Information:

Sequence: Met1-Lys336

Background:

Interleukin3 receptor subunit alpha-2 (IL13RA2/IL3RA2) is also known as also known as cluster of differentiation 213A2 (CD213A2), IL3 receptor subunit alpha-2, IL3R subunit alpha-2, and IL3RA2. The IL13RA2 is often overexpressed in brain tumors, making Il13ra2 one of the vaccine targets for immunotherapy of glioma. IL13RA2/IL3RA2 is a cancer-associated receptor that is present in greater than 80% of High Grade Astrocytomas (HGA) and has recently been recognized as a cytokine that predisposes breast cancer cells to metastasize. Expression of IL13R α 2 was rapidly lost from the surface of transduced cells grown in culture. The loss appeared to be related to ligands present in fetal bovine serum in the medium. None of the malignant glioma cell lines cultivated in vitro and tested to date exhibited the IL13R α 2 receptor. A recombinant virus (R5111) enters cells via its interaction with the IL13R α 2 receptor in a manner that cannot be differentiated from the interaction of wild-type virus with its receptors.