



Recombinant Protein Technical Manual

Recombinant Mouse CD80/B7 Protein (Fc Tag)

RPES2320

Product Data:

Product SKU: RPES2320

Size: 10µg

Species: Mouse

Expression host: Human Cells

Uniprot: Q00609

Protein Information:

Molecular Mass: 50.8 kDa

AP Molecular Mass: 72 kDa

Tag: C-Fc

Bio-activity:

Purity: > 95% as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: T-lymphocyte activation antigen CD80; Activation B7 antigen; B7; CD80;Cd28|Ly-53;Ly53;MIC17;TSA1

Immunogen Information:

Sequence: Val38-Lys245

Background:

Cluster of Differentiation 80, also called B7, is a member of cell surface immunoglobulin superfamily which plays key, yet distinct roles in the activation of T cells. It is the ligand for two different proteins on the T cell surface: CD28 and CTLA-4. Studies have shown that CTLA-4 binds mostly to CD80. The structure presents two extracellular domains: a membrane distal variable-like domain (IgV) and a membrane proximal Ig constant-like domain (IgC) along with an intracellular domain. Both IgV and IgC consist of anti-parallel beta sandwiches joined by a short linker region. CD80 is mostly expressed on the surface of antigen-presenting cells including activated B cells, macrophages and dendritic cells.