



# Recombinant Protein Technical Manual

## Recombinant Human VEGF-D/VEGFD Protein (His Tag)

RPES2504

### Product Data:

**Product SKU:** RPES2504

**Size:** 10µg

**Species:** Human

**Expression host:** Human Cells

**Uniprot:** O43915

### Protein Information:

**Molecular Mass:** 12.2 kDa

**AP Molecular Mass:** 18 kDa

**Tag:** C-6His

**Bio-activity:**

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Endotoxin:** < 1.0 EU per µg as determined by the LAL method.

**Storage:** Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Formulation:** Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

**Reconstitution:** Please refer to the printed manual for detailed information.

**Application:**

**Synonyms:** Vascular Endothelial Growth Factor D; VEGF-D; c-Fos-Induced Growth Factor; FIGF; VEGFD

## Immunogen Information:

**Sequence:** Phe93-Ser201

## Background:

Vascular endothelial growth factor D (VEGF-D) is a member of the platelet-derived growth factor/vascular endothelial growth factor (PDGF/VEGF) family. It is highly expressed in lung, heart, small intestine and fetal lung, and at lower levels in skeletal muscle, colon, and pancreas. VEGF-D is growth factor active in angiogenesis, lymphangiogenesis and endothelial cell growth, stimulating their proliferation and migration and also has effects on the permeability of blood vessels. It may function in the formation of the venous and lymphatic vascular systems during embryogenesis, and also in the maintenance of differentiated lymphatic endothelium in adults. It undergoes a complex proteolytic maturation, generating multiple processed forms that bind and activate VEGFR-2 and VEGFR-3 receptors.