



Recombinant Protein Technical Manual

**Recombinant Human Apolipoprotein D/ApoD
Protein (His Tag)**
RPES2888

Product Data:

Product SKU: RPES2888

Size: 10µg

Species: Human

Expression host: Human Cells

Uniprot: P05090

Protein Information:

Molecular Mass: 20.3 kDa

AP Molecular Mass: 28 kDa

Tag: C-6His

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Apolipoprotein D; Apo-D; ApoD; APOD

Immunogen Information:

Sequence: Glu21-Ser189

Background:

Apolipoprotein-D (ApoD) is an atypical apolipoprotein and, based on its primary structure, it also a member of the lipocalin family. ApoD is mainly associated with high density lipoproteins in human plasma. ApoD is expressed in numerous tissues having high levels of expression in spleen, testes and brain. ApoD plays a role in maintenance and repair within the central and peripheral nervous systems. ApoD occurs in the macromolecular complex with lecithin-cholesterol acyltransferase. It is a multi-ligand, multi-functional transporter and transports a ligand from 1 cell to another. ApoD is probably involved in the transport and binding of bilin, it appears to be able to transport a variety of ligands in a number of different contexts.