



Recombinant Protein Technical Manual
Recombinant Human A1BG/alpha 1B-Glycoprotein
Protein (His Tag)
RPES3274

Product Data:

Product SKU: RPES3274

Size: 50µg

Species: Human

Expression host: HEK293 Cells

Uniprot: NP_570602.2

Protein Information:

Molecular Mass: 53.4 kDa

AP Molecular Mass:

Tag: C-His

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: A1B;ABG;GAB;HYST2477

Immunogen Information:

Sequence: Met 1-Ser495

Background:

AlphaB-glycoprotein (A1BG) is an approximately 80 kDa molecule that circulates in the serum. Mature human A1BG consists of five tandem Ig-like domains and shares 44% amino acid sequence identity with mouse and rat A1BG. Alternative splicing generates an isoform that lacks the signal peptide and the first Ig-like domain. A1BG binds with high affinity to the cysteine-rich protein CRISP-3. In pancreatic carcinoma, it is elevated in pancreatic juice as well as in the serum relative to normal patients. A1BG is abnormally excreted into the urine during idiopathic nephrotic syndrome (INS). Multiple breakdown fragments of A1BG are also present in INS urine, including one of 14 kDa that is found in steroid-resistant but not steroid-sensitive cases. A1BG is also elevated in the urine of patients with interstitial cystitis/painful bladder syndrome (IC/PBS).