



Recombinant Protein Technical Manual

Recombinant Human Fumarase/FH Protein (His Tag)(Active)
RPES3599

Product Data:

Product SKU: RPES3599

Size: 50µg

Species: Human

Expression host: E. coli

Uniprot: P07954

Protein Information:

Molecular Mass: 52 kDa

AP Molecular Mass: 45 kDa

Tag: C-His

Bio-activity: Measured by its ability to transform 1umole of Fumarate to L-malate per minute at pH 7.5 at 37°C. Specific activity is > 25 unit/mg

Purity: > 85 % as determined by reducing SDS-PAGE.

Endotoxin: Please contact us for more information.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile 10mM Tris, 5mM EDTA, 1mM DTT, pH 7.5

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Fumarate Hydratase Mitochondrial; Fumarase; FH;FMRD;HLRCC;LRCC;MCL;MCUL1

Immunogen Information:

Sequence: Ala 45-Lys 510

Background:

Fumarate Hydratase (FH) is an enzymatic component of the tricarboxylic acid (TCA) cycle, or Krebs cycle, and catalyzes the formation of L-malate from fumarate. It exists in both a cytosolic form and an N-terminal extended form, differing only in the translation start site used. The N-terminal extended form is targeted to the mitochondrion, where the removal of the extension generates the same form as in the cytoplasm. Fumarate Hydratase is similar to some thermostable class II fumarases and functions as a homotetramer. Mutations in this gene can cause fumarase deficiency and lead to progressive encephalopathy. Individuals with hemizygous germline fumarate hydratase (FH) mutations are predisposed to renal cancer. These tumors predominantly exhibit functional inactivation of the remaining wild-type allele, implicating FH inactivation as a tumor-promoting event.