



# Recombinant Protein Technical Manual

## Recombinant Rat CD40/TNFRSF5 Protein (His Tag)

RPE5108

### Product Data:

**Product SKU:** RPE5108

**Size:** 50µg

**Species:** Rat

**Expression host:** HEK293 Cells

**Uniprot:** NP\_599187.1

### Protein Information:

**Molecular Mass:** 20.6 kDa

**AP Molecular Mass:** 28 kDa

**Tag:** C-His

**Bio-activity:**

**Purity:** > 95 % as determined by SDS-PAGE

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method

**Storage:** Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Formulation:** Lyophilized from sterile PBS, pH 7.4

**Reconstitution:** Please refer to the printed manual for detailed information.

**Application:**

**Synonyms:** CD40;Tnfrsf5

## Immunogen Information:

**Sequence:** Met1-Arg193

## Background:

CD40, also known as TNFRSF5, is a member of the TNF receptor superfamily which are single transmembrane-spanning glycoproteins. CD40 protein plays an essential role in mediating a broad variety of immune and inflammatory responses including T cell-dependent immunoglobulin class switching, memory B cell development, and germinal center formation. CD40 protein is expressed in B cells, dendritic cells, macrophages, endothelial cells, and several tumor cell lines. Defects in CD40 result in hyper-IgM immunodeficiency type 3 (HIGM3). In addition, CD40/CD40L interaction is found to be necessary for amyloid-beta-induced microglial activation, and thus is thought to be an early event in Alzheimer disease pathogenesis.