

Mouse IgG2a, kappa Isotype Control [C1.18.4] -Low Endotoxin Purified AGEL2582

Description

This Mouse IgG2a, kappa Isotype Control [C1.18.4] -Low Endotoxin Purified is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL2582
Contents:	100µg, 25µg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	C1.18.4
Synonyms:	-
Applications:	FCM
Reactivity:	-
Immunogen:	-

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	-
Host Species:	Mouse
Isotype:	Mouse IgG2a, κ
Isotype Control:	-

Conjugation:	-
Conjugation Information:	-
Buffer:	Sterile PBS, pH 7.2. < 1.0 EU per mg of the antibody as determined by the LAL method.
Purification:	-
Target:	IgG2a
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	-
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles. This preparation contains no preservatives, thus it should be handled under aseptic conditions.

Shipping: Ice bag

Recommended Dilution: -

Recommended Usage:

Application	Recommended Usage
-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.