

Purified Anti-Mouse/Human GAL3 Antibody [M3/38]

AGEL5137

Description

This Purified Anti-Mouse/Human GAL3 Antibody [M3/38] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL5137
Contents:	25µg, 100µg, 1mg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	M3/38
Synonyms:	MAC, Galactose-specific lectin, Galactoside-BindingSoluble, Galectin, Lectin L, GAL, CBP, LGALS3, CBP35, GAL3, GALBP, GALIG, L31, LGALS2, MAC2, Galactoside-BindingSoluble 3, Lectin, Galectin-3, Gal-3, L-31, 35 kDa lectin, CBP 35, Galactose-specific lectin 3, IgE-binding protein, Laminin-binding protein, Lectin L-29, Mac-2 antigen, Carbohydrate binding protein 35, Carbohydrate-binding protein 35, galactose binding, galactoside bindingSoluble 3, Galactoside-binding protein, galectin 3, Galectin 3 internal gene, Galectin3, GBP, IgE binding protein, included, L 31, L 34, L-34 galactoside-binding lectin, Lectin, Galactoside-BindingSoluble 3, LEG3, MAC 2 antigen, Mac-2, Macrophage galactose-specific lectin, MGC105387, P35,Soluble 3
Applications:	FCM
Reactivity:	Mouse/Human
Immunogen:	Recombinant Mouse GAL3 protein

Antibody Data

Uniprot ID: -

Gene ID:	-
Swissprot:	P17931
Host Species:	Rat
Isotype:	Rat IgG2a, κ
Isotype Control:	-
Conjugation:	-
Conjugation Information:	-
Buffer:	Phosphate-buffered solution, pH 7.2, containing 0.05% non-protein stabilizer. Dialyze to completely remove the stabilizer prior to labeling.
Purification:	>98%, Protein A/G purified
Target:	GAL3
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	Verified Samples in FCM: Human peripheral blood lymphocytes
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles.

Shipping: Ice bag

Recommended Dilution: FCM 2 µg/mL (0.5×10⁶-1×10⁶ cells)

Recommended Usage:

Application	Recommended Usage
-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes:

Centrifuge before opening to ensure complete recovery of vial contents.