

Anti-RANK [R07-1T2] Monoclonal Antibody

AGMB00017

Description

This Anti-RANK [R07-1T2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00017
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD265, FEO, LOH18CR1, mRANK, ODFR, OFE, OPTB7, OSTs, PDB2, RANK, TNFRSF11A, TRANCER
Applications:	
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-1T2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9Y6Q6
Immunogen:	A synthesized peptide derived from human RANK
Gene ID:	8792

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: TNFRSF11A
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 66 kDa, Observed MW: 80 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.