

Anti-MRPS26 [R01-2D6] Monoclonal Antibody

AGMB00299

Description

This Anti-MRPS26 [R01-2D6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00299
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	C20orf193, GI008, MRP S13, MRP S26, MRPS13, MRPS26, NY BR 87, RPMS13, S13mt, S26mt
Applications:	WB ICC/IF
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R01-2D6
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9BYN8
Immunogen:	Recombinant protein of human MRPS26
Gene ID:	64949

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: MRPS26
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 24 kDa, Observed MW: 24 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	ICC/IF	1:100-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	