

Anti-CUG BP1 [R03-0G7] Monoclonal Antibody

AGMB00452

Description

This Anti-CUG BP1 [R03-0G7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00452
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CELF1, BRUNOL2, CUGBP, CUGBP1, NAB50, CUGBP Elav-like family member 1, CELF-1, 50 kDa nuclear polyadenylated RNA-binding protein, Bruno-like protein 2, CUG triplet repeat RNA-binding protein 1, CUG-BP1, CUG-BP- and ETR-3-like factor 1, Dead, 50 kDa Nuclear polyadenylated RNA binding protein, 50 kDa nuclear polyadenylated RNA-binding protein, Bruno like 2,
Applications:	WB IHC-P ICC/IF
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R03-0G7
Clonality:	Monoclonal Antibody
SwissProt ID:	Q92879

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Immunogen: A synthetic peptide of human CUG-BP1

Gene ID: 10658

Gene Name: CELF1

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 52 kDa, Observed MW: 52 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC	1:50-1:100
	ICC/IF	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.