

Anti-PPP2R5E [R06-6E2] Monoclonal Antibody

AGMB00677

Description

This Anti-PPP2R5E [R06-6E2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00677
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	B56E, B56epsilon, 2A5E_HUMAN, Epsilon isoform of regulatory subunit B56 protein phosphatase 2A, PP2A B subunit B' epsilon,
Applications:	WB IHC-F IHC-P ICC/IF FC IP
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-6E2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q16537
Immunogen:	A synthetic peptide of human PPP2R5E
Gene ID:	5529

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	PPP2R5E
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 55 kDa, Observed MW: 55 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:200-1:500
	ICC/IF	1:200-1:500
	FC	1:200-1:1000
	IP	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	