

Anti-LLGL1 [R05-0B8] Monoclonal Antibody

AGMB00722

Description

This Anti-LLGL1 [R05-0B8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00722
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	DLG4, HUGL, LLGL, Lgl1, Mgl1, HUGL1, HUGL-1
Applications:	WB IHC-P ICC/IF IP
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-0B8
Clonality:	Monoclonal Antibody
SwissProt ID:	Q15334
Immunogen:	Recombinant protein of human LLGL1
Gene ID:	3996
Gene Name:	LLGL1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 115 kDa, Observed MW: 115 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:50-1:100
	ICC/IF	1:50-1:100
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.