

Anti-Flotillin 2 [R01-3I1] Monoclonal Antibody

AGMB00784

Description

This Anti-Flotillin 2 [R01-3I1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00784
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Epidermal Surface Antigen, ESA, Flotillin-2, Membrane Component Chromosome 17 Surface Marker 1, ESA1, M17S1, FLOT2, 1E11G12, ECS 1
Applications:	
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-3I1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q14254
Immunogen:	A synthesized peptide derived from human Flotillin 2
Gene ID:	2319

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: FLOT2
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 47 kDa, Observed MW: 47 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IP	1:30
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	