

Anti-HMGCS2 [R04-2C5] Monoclonal Antibody

AGMB01221

Description

This Anti-HMGCS2 [R04-2C5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01221
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	HMG CoA synthase, HMGCS2, hydroxy 3 methylglutaryl Coenzyme A synthase 2
Applications:	WB IHC-P
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R04-2C5
Clonality:	Monoclonal Antibody
SwissProt ID:	P54868
Immunogen:	A synthesized peptide derived from human HMGCS2
Gene ID:	3158

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	HMGCS2
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 57 kDa, Observed MW: 50 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	