

Anti-Cleaved-Caspase 2 [R04-3F4] Monoclonal Antibody

AGMB01421

Description

This Anti-Cleaved-Caspase 2 [R04-3F4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01421
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Caspase-2, Caspase 2, Caspase2, CASP-2, Neural precursor cell expressed developmentally down-regulated protein 2, NEDD-2, Protease ICH-1, CASP2, ICH1, NEDD2
Applications:	WB IHC-P ICC/IF FC
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R04-3F4
Clonality:	Monoclonal Antibody
SwissProt ID:	P42575
Immunogen:	A synthesized peptide derived from human Caspase-2
Gene ID:	835

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	CASP2
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 51 kDa, Observed MW: 48,30,12 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.