

Anti-CCR2 [R03-5F4] Monoclonal Antibody

AGMB01439

Description

This Anti-CCR2 [R03-5F4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01439
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CC-CKR-2, CCR-2, CCR2A, CCR2B, CD192, CKR2, CKR2A, CKR2B, CMKBR2, MCP-1-R
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	R03-5F4
Clonality:	Monoclonal Antibody
SwissProt ID:	P41597
Immunogen:	A synthesized peptide derived from human CCR2/CKR2
Gene ID:	729230

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CCR2
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 42kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:100-1:200
	IP	1:50
	FC	1:50-1:100
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.