

Anti-Aryl hydrocarbon Receptor [R07-1H9] Monoclonal Antibody

AGMB01497

Description

This Anti-Aryl hydrocarbon Receptor [R07-1H9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01497
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	AhR, Ah receptor, Class E basic helix-loop-helix protein 76, bHLHe76, AHR, RP85, bHLHe76, AHR
Applications:	WB ICC/IF
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-1H9
Clonality:	Monoclonal Antibody
SwissProt ID:	P35869
Immunogen:	A synthetic peptide of human Aryl hydrocarbon Receptor
Gene ID:	196

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	AHR
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 96 kDa, Observed MW: 100 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	ICC/IF	1:100-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	