

Anti-GCK [R04-7V9] Monoclonal Antibody

AGMB01512

Description

This Anti-GCK [R04-7V9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01512
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	GCK, Glucokinase, Hexokinase type IV, HK IV, Hexokinase-4, HK4, Hexokinase-D
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R04-7V9
Clonality:	Monoclonal Antibody
SwissProt ID:	P35557
Immunogen:	A synthesized peptide derived from human Glucokinase
Gene ID:	2645

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	GCK
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 52 kDa, Observed MW: 52 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.