

Anti-DNAJA1 [R02-7J8] Monoclonal Antibody

AGMB01589

Description

This Anti-DNAJA1 [R02-7J8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01589
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	DNAJA1, DNAJ2, HDJ2, HSJ2, HSPF4, DnaJ homolog subfamily A member 1, DnaJ protein homolog 2, HSDJ, Heat shock 40 kDa protein 4, Heat shock protein J2, HSJ-2, Human DnaJ protein 2, hDj-2
Applications:	WB FC IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R02-7J8
Clonality:	Monoclonal Antibody
SwissProt ID:	P31689
Immunogen:	A synthesized peptide derived from human DNAJA1
Gene ID:	3301

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	DNAJA1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 45 kDa, Observed MW: 45 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	FC	1:50-1:100
	IP	1:50
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	