

# Anti-Cyclin D3/CCND3 [R07-6K2] Monoclonal Antibody

## AGMB01591

### Description

---

This Anti-Cyclin D3/CCND3 [R07-6K2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                  |
|-----------------------|------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB01591                                                        |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)              |
| <b>Synonyms:</b>      | CCND 3, Ccnd3, CyclinD3, D3 type cyclin, G1 S specific cyclin D3 |
| <b>Applications:</b>  | <b>WB</b> <b>ICC/IF</b> <b>FC</b>                                |
| <b>Research Area:</b> | Cell Biology                                                     |
| <b>Form:</b>          | Liquid                                                           |

### Antibody Data

---

|                      |                                                    |
|----------------------|----------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse                                       |
| <b>Clone:</b>        | R07-6K2                                            |
| <b>Clonality:</b>    | Monoclonal Antibody                                |
| <b>SwissProt ID:</b> | P30281                                             |
| <b>Immunogen:</b>    | A synthesized peptide derived from human Cyclin D3 |
| <b>Gene ID:</b>      | 896                                                |
| <b>Gene Name:</b>    | CCND3                                              |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Host Species:</b>     | Rabbit                                     |
| <b>Isotype:</b>          | IgG                                        |
| <b>Purification:</b>     | Affinity Chromatography                    |
| <b>Conjugated:</b>       | Unconjugated                               |
| <b>Modification:</b>     | Unmodified                                 |
| <b>Molecular Weight:</b> | Calculated MW: 33 kDa, Observed MW: 33 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| <b>Antibody Dilution Ratio:</b> | Application | Antibody Dilution Ratio |
|---------------------------------|-------------|-------------------------|
|                                 | WB          | 1:1000-1:2000           |
|                                 | IF          | 1:50-1:200              |
|                                 | FC          | 1:50-1:100              |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.