

Anti-MST1 [R01-4B7] Monoclonal Antibody

AGMB01657

Description

This Anti-MST1 [R01-4B7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01657
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MSP, HGFL, NF15S2, D3F15S2, DNF15S2, MSP/MST1
Applications:	WB IHC-P FC
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-4B7
Clonality:	Monoclonal Antibody
SwissProt ID:	P26927
Immunogen:	A synthesized peptide derived from human MST1
Gene ID:	4485
Gene Name:	MST1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 80 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.