

Anti-Galectin 3 [R05-8J4] Monoclonal Antibody

AGMB01696

Description

This Anti-Galectin 3 [R05-8J4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01696
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LGALS3, 35 kDa lectin, Carbohydrate binding protein 35, Carbohydrate-binding protein 35, Galectin-3, Gal-3, CBP 35l, GALBP, Galactose-specific lectin 3, Galactoside-binding protein (GALBP)
Applications:	WB IHC-P ICC/IF FC
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R05-8J4
Clonality:	Monoclonal Antibody
SwissProt ID:	P17931
Immunogen:	A synthetic peptide of human Galectin 3
Gene ID:	3958

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: LGALS3

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 26 kDa, Observed MW: 26 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:100-1:200
	ICC/IF	1:500-1:1000
	FC	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.