

Anti-C1QA [R29-1P-7] Monoclonal Antibody

AGMB01794

Description

This Anti-C1QA [R29-1P-7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01794
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	C1qa, Complement C1q subcomponent subunit A, Complement component 1 q subcomponent A chain, Complement component 1 q subcomponent alpha polypeptide, C1QA_HUMAN, C1QB,
Applications:	WB IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human, Rat
Clone:	R29-1P-7
Clonality:	Monoclonal Antibody
SwissProt ID:	P02745
Immunogen:	A synthesized peptide derived from human C1QA
Gene ID:	712

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: C1QA
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW:26kDa, Observed MW:30kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:1000-1:2000
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.