

Anti-HTF9C/TRMT2A [R03-3H4] Monoclonal Antibody

AGMB01938

Description

This Anti-HTF9C/TRMT2A [R03-3H4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01938
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	HpaII tiny fragments locus 9c protein, HpaII tiny fragments locus 9c protein, HTF9C, TRM2 tRNA methyltransferase 2 homolog A (S. cerevisiae),
Applications:	WB IHC-F IHC-P ICC/IF
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-3H4
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8IZ69
Immunogen:	A synthetic peptide of human HTF9C
Gene ID:	27037

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	TRMT2A
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 69 kDa, Observed MW: 69 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:50-1:100
	ICC/IF	1:100-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	