

Anti-Heme Oxygenase 1 [R03-8A2] Monoclonal Antibody

AGMB01960

Description

This Anti-Heme Oxygenase 1 [R03-8A2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01960
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	HO-1, HSP32, HMOX1D, bK286B10, Heme Oxygenase 1 (HO-1/HMOX1)
Applications:	WB ICC/IF FC IP
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R03-8A2
Clonality:	Monoclonal Antibody
SwissProt ID:	P09601
Immunogen:	A synthetic peptide of human Heme Oxygenase 1
Gene ID:	3162
Gene Name:	HMOX1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 33 kDa, Observed MW: 33 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	ICC/IF	1:100-1:200
	FC	1:100-1:200
	IP	1:20

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.