

Anti-TRIM23/ARD1 [R09-4K7] Monoclonal Antibody

AGMB02388

Description

This Anti-TRIM23/ARD1 [R09-4K7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02388
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ARD1, ARFD1, RNF46, ADP ribosylation factor domain protein 1, ADP-ribosylation factor domain-containing protein 1, ARF domain protein 1,
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R09-4K7
Clonality:	Monoclonal Antibody
SwissProt ID:	P36406
Immunogen:	A synthetic peptide of human TRIM23
Gene ID:	373

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	TRIM23
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 64 kDa, Observed MW: 64 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:2000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.