

Anti-TXNL1 [R06-8C2] Monoclonal Antibody

AGMB02393

Description

This Anti-TXNL1 [R06-8C2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02393
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	TRP32, TXL, Thioredoxin-like protein 1, Thioredoxin like protein 1, thioredoxin like 1
Applications:	WB IHC-P ICC/IF FC
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-8C2
Clonality:	Monoclonal Antibody
SwissProt ID:	O43396
Immunogen:	A synthetic peptide of human TXNL1
Gene ID:	9352

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	TXNL1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 32 kDa, Observed MW: 32 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	