

Anti-ZBTB7A [R07-5E0] Monoclonal Antibody

AGMB02569

Description

This Anti-ZBTB7A [R07-5E0] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02569
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FBI1, LRF, Pokemon, TIP21, ZBTB7, ZBTB7A, ZNF857A
Applications:	WB IHC-P ICC/IF IP
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-5E0
Clonality:	Monoclonal Antibody
SwissProt ID:	O95365
Immunogen:	A synthetic peptide of human ZBTB7A
Gene ID:	51341
Gene Name:	ZBTB7A

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 61 kDa, Observed MW: 61 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:100-1:200
	ICC/IF	1:200-1:500
	IP	1:10-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	