

# Anti-PKC beta 2 [R06-0C2] Monoclonal Antibody

## AGMB02571

### Description

---

This Anti-PKC beta 2 [R06-0C2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                  |
|-----------------------|------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB02571                                                        |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)              |
| <b>Synonyms:</b>      | PRKCB, PKCB, PRKCB1, Protein kinase C beta type, PKC-B, PKC-beta |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b> <b>ICC/IF</b> <b>FC</b>                   |
| <b>Research Area:</b> | Signal Transduction                                              |
| <b>Form:</b>          | Liquid                                                           |

### Antibody Data

---

|                      |                                                     |
|----------------------|-----------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                                   |
| <b>Clone:</b>        | R06-0C2                                             |
| <b>Clonality:</b>    | Monoclonal Antibody                                 |
| <b>SwissProt ID:</b> | P05771                                              |
| <b>Immunogen:</b>    | A synthesized peptide derived from human PKC beta 2 |
| <b>Gene ID:</b>      | 5579                                                |
| <b>Gene Name:</b>    | PRKCB                                               |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Host Species:</b>     | Rabbit                                     |
| <b>Isotype:</b>          | IgG                                        |
| <b>Purification:</b>     | Affinity Purified                          |
| <b>Conjugated:</b>       | Unconjugated                               |
| <b>Modification:</b>     | Unmodified                                 |
| <b>Molecular Weight:</b> | Calculated MW: 77 kDa, Observed MW: 77 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | WB                 | 1:500-1:1000                   |
|                                 | IHC-P              | 1:200-1:500                    |
|                                 | ICC/IF             | 1:50-1:100                     |
|                                 | FC                 | 1:20-1:50                      |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.