

Anti-Ly-6A/E (Sca-1) [R05-8A2] Monoclonal Antibody

AGMB02703

Description

This Anti-Ly-6A/E (Sca-1) [R05-8A2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02703
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	TAP, Sca1, Sca-1, Ly-6A.2, Ly-6A/E, Ly-6E.1, Ly-6A/E (Sca-1)
Applications:	WB IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R05-8A2
Clonality:	Monoclonal Antibody
SwissProt ID:	P05533
Immunogen:	A synthesized peptide derived from mouse LY6A
Gene ID:	110454
Gene Name:	Ly6A

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 14 kDa, Observed MW: 10-20 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.