

Anti-NEAS [R01-9G3] Monoclonal Antibody

AGMB02899

Description

This Anti-NEAS [R01-9G3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02899
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	SPTAN1, NEAS, SPTA2, Spectrin alpha chain, non-erythrocytic 1, Alpha-II spectrin, Fodrin alpha chain, Spectrin, non-erythroid alpha subunit, (ALPHA)II-SPECTRIN, Alpha II Spectrin, Alpha-II spectrin,
Applications:	WB IHC-P IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-9G3
Clonality:	Monoclonal Antibody
SwissProt ID:	Q13813
Immunogen:	A synthesized peptide derived from human Alpha Fodrin
Gene ID:	6709

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: SPTAN1

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 285 kDa, Observed MW: 150,285 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IP	1:10

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.