

Anti-GSTM1 [R05-1W9] Monoclonal Antibody

AGMB03145

Description

This Anti-GSTM1 [R05-1W9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03145
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Glutathione S-transferase Mu 1, GSTM1, MU, H-B, GST1, GTH4, GTM1, MU-1, GSTM1-1, MGC26563, GSTM1a-1a, GSTM1b-1b
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R05-1W9
Clonality:	Monoclonal Antibody
SwissProt ID:	P09488
Immunogen:	A synthesized peptide derived from human Glutathione S Transferase mu
Gene ID:	2944

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: GSTM1

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 26 kDa, Observed MW: 26 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.