

Anti-HOXA13 [R01-1S7] Monoclonal Antibody

AGMB03195

Description

This Anti-HOXA13 [R01-1S7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03195
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Homeobox1J, HomeoboxA13, HOX 1, HOX 1J, HOX A13, HOX1, HOX1J, HOXA 13, HOXA13, Transcription factor HOXA13
Applications:	
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-1S7
Clonality:	Monoclonal Antibody
SwissProt ID:	P31271
Immunogen:	A synthesized peptide derived from human HOXA13
Gene ID:	3209

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	HOXA13
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 40 kDa, Observed MW: 40 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.