

Anti-HLA DMB [R24-3U-2] Monoclonal Antibody

AGMB03396

Description

This Anti-HLA DMB [R24-3U-2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03396
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	class II antigen HLA-DM beta chain, class II histocompatibility antigen, Class II histocompatibility antigen M beta chain,
Applications:	
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R24-3U-2
Clonality:	Monoclonal Antibody
SwissProt ID:	P28068
Immunogen:	A synthesized peptide derived from human HLA DMB
Gene ID:	3109

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: HLA-DMB

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 29 kDa, Observed MW: 29 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	FC	1:50
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	