

Anti-ULBP1 [R86-8D-9] Monoclonal Antibody

AGMB03452

Description

This Anti-ULBP1 [R86-8D-9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|-----------------------|--|
| SKU: | AGMB03452 |
| Contents: | 20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml) |
| Synonyms: | Alcan beta, N2DL1, NKG2DL1, RAET1I, UL16, UL16 binding protein 1, ULBP1, |
| Applications: | WB |
| Research Area: | Immunology |
| Form: | Liquid |

Antibody Data

| | |
|----------------------|--|
| Reactivity: | Human |
| Clone: | R86-8D-9 |
| Clonality: | Monoclonal Antibody |
| SwissProt ID: | Q9BZM6 |
| Immunogen: | A synthesized peptide derived from human ULBP1 |
| Gene ID: | 80329 |

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ULBP1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW:28kDa, Observed MW:36kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|---------------------------------|--------------------|--------------------------------|
| | WB | 1:500-1:1000 |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.