

Anti-Cathepsin H [R53-5Z-4] Monoclonal Antibody

AGMB03460

Description

This Anti-Cathepsin H [R53-5Z-4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03460
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ACC4, ACC5, CPSB, ACC-4, ACC-5, CTSH, CPSB, Pro-cathepsin H
Applications:	WB IHC-P ICC/IF
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R53-5Z-4
Clonality:	Monoclonal Antibody
SwissProt ID:	P09668
Immunogen:	A synthesized peptide derived from human Cathepsin H
Gene ID:	1512
Gene Name:	CTSH

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW:37kDa, Observed MW:37kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:100-1:200
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	