

# Anti-Cyclin D3/CCND3 [R05-9X-5] Monoclonal Antibody

## AGMB03587

### Description

---

This Anti-Cyclin D3/CCND3 [R05-9X-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB03587
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	CCND 3, Ccnd3, CyclinD3, D3 type cyclin, G1 S specific cyclin D3
<b>Applications:</b>	<b>IHC-P</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R05-9X-5
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P30281
<b>Immunogen:</b>	A synthesized peptide derived from human Cyclin D3
<b>Gene ID:</b>	896
<b>Gene Name:</b>	CCND3

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	IHC	1:100-1:200

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.