

Anti-SSRP1 [R01-7Z-1] Monoclonal Antibody

AGMB03915

Description

This Anti-SSRP1 [R01-7Z-1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03915
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FACT, T160, FACT80, SSRP1, Chromatin-specific transcription elongation factor 80 kDa subunit, Cisplatin DNA SSRP, Facilitates chromatin remodeling 80 kDa subunit,
Applications:	WB
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-7Z-1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q08945
Immunogen:	Peptide
Gene ID:	6749

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: SSRP1

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW:81 kDa, Observed MW: 81 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.