

Anti-BTN3A1 [R05-7Z-5] Monoclonal Antibody

AGMB03957

Description

This Anti-BTN3A1 [R05-7Z-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03957
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	BTN3A1, BTF5, CD277, BTN3.1, BT3.1, BT3A1_HUMAN, Butyrophilin subfamily 3 member A1
Applications:	FC ELISA
Research Area:	Immunology
Form:	Lyophilized

Antibody Data

Reactivity:	Human
Clone:	R05-7Z-5
Clonality:	Monoclonal Antibody
SwissProt ID:	O00481
Immunogen:	Recombinant human BTN3A1 protein
Gene ID:	11119

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: BTN3A1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at -20°C to -80°C for 12 months in lyophilized form.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	FC	1:50-1:100
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.