

Anti-ERG [R09-9J7] Monoclonal Antibody

AGMB04069

Description

This Anti-ERG [R09-9J7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04069
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ERG, Transcriptional regulator ERG, Transforming protein ERG, Avian erythroblastosis virus E-26 (v-ets) oncogene related, D030036I24Rik, Erg 3,
Applications:	WB IHC-P ICC/IF FC
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R09-9J7
Clonality:	Monoclonal Antibody
SwissProt ID:	P11308
Immunogen:	A synthesized peptide derived from human ERG
Gene ID:	2078

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ERG
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 54 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:200-1:1000
	ICC/IF	1:200-1:1000
	FC	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.